Basic Concepts and Materials used in Unilateral Biportal Endoscopy

Gonçalo Januário

Neurosurgeon Consultant, Neurosurgery Department, Policlínica Juaneda Miramar, Palma de Mallorca, Balearic Islands, Spain.

*Corresponding Author: Gonçalo Januário, Neurosurgeon Consultant, Neurosurgery Department, Policlínica Juaneda Miramar, Palma de Mallorca, Balearic Islands, Spain.

Received date: January 09, 2025; Accepted date: January 24, 2025; Published date: February 05, 2025

Citation: Gonçalo Januário, (2025), Basic Concepts and Materials used in Unilateral Biportal Endoscopy, *J. Neuroscience and Neurological Surgery*, 17(2); **DOI:10.31579/2578-8868/358**

Copyrights: ©, 2025, Gonçalo Januário. This is an open-access article distributed under the terms of The Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Introduction

The French physician Antoine Jean Desormeaux, lived in the 19thcentury, considered the "father of endoscopy" developed an instrument to examine the urinary tract and bladder named "endoscope". It was the first time in history that this term was used.1

Throughout the last century, many technological advances were achieved with increasingly sophisticated devices and surgical approaches. Particularly during the last decade the endoscopic surgery began to increase the applicability for the treatment of spine pathologies.

Actually is not yet considered gold standard but is progressing quickly to become the reference technique maybe during the next years. Like all surgical techniques present advantages and disadvantages, ideal for decompression and discectomy in cases of stenosis or disc herniation respectively. At this time it is increasingly used in the 3 segments of the spine (cervical, dorsal, lumbar) and in different pathologies. However, we must not forget that it has a limited indication in cases of high-grade listhesis as an isolated technique.

Two endoscopic surgical techniques are currently available, uniportal and biportal. Uniportal, full- endoscopic spine surgery is a minimally invasive procedure that uses a keyhole approach for the decompression or the discectomy. For other part the Unilateral biportal endoscopy (UBE), can be performed more widely, use two portals that permit high range of movements and amplified vision in the surgical field. When compared both techniques the muscle destruction is similar.2

All the endoscopic procedures require a working space, but the spine from an anatomical point of view present a limited work space especially by the bone and osteoarticular structures.

A important meaning in the endoscopic spine surgeries is the Kambin's triangle, an anatomical reference, a natural corridor that allows the ideal approach to the lumbar canal.3

As in all spine surgeries, also in both endoscopic techniques, the manipulation of the duramater and particularly the nerve roots should be avoided to reduce complications.

The UBE requires two or more channels, in cases with a multilevel approach, to create sufficient work spaces around the levels we want to treat. The approaches in UBE are interlaminar, between the two laminae, or extraforaminal, depending the pathology location.

The anatomical corridor for the interlaminar approach, present two potential ways. The interfascicular related with the two small muscles of the multifidus, and the space containing fat and connective tissue located between the multifidus and the lamina. When these two potential spaces can be converted to an atraumatic working path, UBE guarantees the same indications as conventional spinal surgical treatments but with concepts of minimally invasive surgery.

In our center we don't use dilators to do the approach and find the interlaminar space, we respect the anatomical corridors. For the extraforaminal approach, the surrounding area of the foramen with fat and connective tissue acts as the working space.4

As a big difference for open surgery the endoscopic spinal surgery is a fluid-medium procedure, similar to orthopedic arthroscopy. The operational concept is similar in both procedures.

In UBE we need to create two surgical portals on the same side of the spine, one for the endoscope and the other for the instruments used during the procedure.

The irrigating with saline solution enters in the endoscopic portal and leaves through the working portal (Figure 1).




Figure 1: illustration of fluid input and output in UBE.

A continuous and controlled circulation of the saline solution fluid is a critical factor for performing a successful UBE. The ideal pressure in the working space is 30–50 mmHg. The correct circulation of the saline solution between the channels allows a clearer view of the surgical field and in addition reduce the infection rate.

Therefore, to prevent complications related to circulation and working in a liquid medium, which are one of the main problems for beginners, it's important facilitate the free circulation of the flow between the portals. In UBE several structures can influence and limit the free flow of saline solution, such as the fascia and the posterior paravertebral muscles. Furthermore, the spinal canal is surrounded by bony and ligamentous structures that also reduce the free circulation of the flow. Once the yellow ligament is removed, in the interlaminar approaches, the epidural space is directly entered without a separating structure in addition to the thin layer of epidural fat. Attention in this point for evitate dural tears or leaks As in all surgical techniques, knowledge of anatomy and anatomical landmarks is essential. It allows a better understanding and guidance during the procedure, avoiding complications and facilitating the workflow. It certainly reduces the need for intraoperative x-rays or the use of the neuronavigation system. Mastery of anatomical references can allow the surgeon to realize that the navigation system may be presenting a shift on relation to the real location. All the technology available in the operating room is a great advantage but the degree of accuracy must always be verified at all times. This is an important point for all those who begin to perform this technique.

Another key point that we must take into account is the distance from the skin to the spine. This aspect, especially in obese or with highly developed muscles patients, makes it difficult for the flow of saline solution to

circulate. And it can be a confusion factor or complications during surgery, especially in the first steps of the surgical procedure. If this happens, do not forget the basic anatomical points, identifying them at all times and of course confirming them with the intraoperative imaging techniques that we have available in our operating rooms. Do not forget that it is very important confirm the correct opening of the fascia.

If the circulation is not correct or the pressure increases a high hydrostatic pressure may cause headache in the postoperative period. To prevent complications the use of a controlled pressure pump or using the effect of gravity with the principle of communicating vessels is recommended. Thus avoiding high intracanal pressure in our work space. The optimal hydrostatic pressure, depend the authors, is located between 30–50 Auctores Publishing LLC – Volume 17(2)-358 www.auctoresonline.org ISSN: 2578-8868

mmHg, which varies according to the distance between the irrigating saline bag and the working space. The working space in the patient is fixed less then 1 cm and the height of the saline bag decides the pressure.

For calculate the pressure (mmHg) we dividing the distance (cm) between the patient and the saline bag by 1.36. [height(cm) \div 1.36 = pressure (mmHg)].[5]

In our center we don't use pump-controlled system, we use the effect of the gravity for a natural input and output of the saline solution. Sometimes, if necessary, we temporarily place a drainage to help the output flow. One important point for facilitate the flow between the portals, in our opinion, is a correct open of the fascia.

The principal factors related to fluid complications are a excessive saline pressure in the epidural space impacts the neural system and can cause headache, muscle pain, transient dysesthesia or paresthesias, nuchal pain, long-level epidural bleeding, delayed awakening from general anesthesia, seizure, intracranial hematoma, blindness due to retinal bleeding.

Are also described as a rare complication like a massive abdominal fluid collection could occur while performing the paraspinal approach. Is caused by massive damage to the inter-transverse membrane and ligaments, which function as a barrier between the back and abdominal cavity.

Another complication related with the fluid could be hypothermia caused by the prolonged use of cold saline, for prevention we recommend a tepid saline and a heating pad on the surgical table.

Materials in UBE:

1. Endoscope:

- High-Definition Video Endoscope: Provides clear and magnified visualization of the surgical site.
- Rigid or Flexible Endoscope: Typically, rigid endoscopes are used for spine surgeries, as they provide stability and precise visualization.

2. Surgical Instruments:

- Specialized Endoscopic Instruments: Includes graspers, scissors, dissectors, and burrs designed to fit through small incisions.
- Motorized Tools: Some UBE procedures may involve the use of motorized devices, such as drills or shavers, to assist with tissue

J. Neuroscience and Neurological Surgery

removal or bone resection.

3. Trocars and Cannulas: in our center we don't use, but we make reference.

- Cannulas: Hollow tubes inserted into the body to allow for instrument passage and to maintain a working space.
- Trocars: Sharp devices used to create the initial access point for the cannulas.

4. Surgical Navigation and Imaging:

- Intraoperative Navigation Systems: Provide real-time guidance to ensure accurate placement of instruments and reduce risks of nerve damage.
- Fluoroscopy (X-ray Imaging): Used in some cases for real-time imaging, especially during spinal procedures to check for alignment and instrumentation placement.

5. Coagulation devices:

- Electrocautery: Used to cauterize tissue and control bleeding.
- Bipolar coagulation device.
- Radiofrequency system, unavailable in our center.
- Lasers: Sometimes used for tissue removal or vaporization, depending on the procedure, unavailable in our center. Unavailable in our center.

6. Materials for decompression, discectomy:

- Kerrison punch.
- Curettes.

- Dissectors
- Indian Knife
- Used to keep the surgical area clean, flush out debris, and maintain a clear view during the procedure.

7. Saline Solution or Irrigation System:

• Used to keep the surgical area clean, flush out debris, and maintain a clear view during the procedure.

8. Suture Materials:

• Endoscopic Suturing Tools: If closure is needed for wounds or tissue, special suturing tools designed for minimally invasive techniques are used.

The UBE approach can be done in a operation room with a team composed for surgeon and a assistant, anesthesiologist, scrub nurse and circulating nurse. Commonly we use the endoscope of 0° with the benefits that this lens reduces the vision distortion that could be caused by an angled lens. Sometimes we need to use also the 30° degrees angulation view especially in case of bilateral decompression and to observe the dead angles. The other characteristics of the endoscope are 176 mm length and 4 mm diameter. During the approach the drilling system is very helpful, we use the drill burr with a sheath as in arthroscopy. However, occasionally, a high-speed drill without a sheath is applied for areas near sensitive cord levels. In our center we use a standard endoscopic drill burr with a sheath, could be effective for controlling the water pressure in the working space and for bone dust drainage.

For the procedure we use as equipment: endoscopic system (endoscope, camera, cable, and monitor) and drilling system (Figure 2).




Figure 2: endoscopic (A,B,C) and drilling system (D,E).

For localize the level to operate we use a C-arm fluoroscopy system or neuronavigation, we use the fusion of both system as a guide to select the correct entry points and also to accurate the position/location during surgery (Figure 3).

J. Neuroscience and Neurological Surgery

Copy rights @ Gonçalo Januário,



Figure 3: A: c-arm device; B: neuronavigation system; C: intraoperative image with fusion or radiography and neuronavigation.

In UBE surgery we use two different types of coagulation systems monopolar and bipolar. The monopolar coagulation is used in all the procedures except in patient with a pacemaker.

The irrigation type is also very important, the saline solution is a very effective source and cooperate in the coagulation process with greater conductivity.

For control the hemostasis and maintain a clear vision we use monopolar device outside the medular canal, inside the canal and near the duramater or nerve roots we use bipolar coagulation.

Outside the canal, during the first steps of the procedure, with a intensity between 15-25, the monopolar system (for muscles, bone and soft tissues)

permit a good bleeding control. Inside the canal we use bipolar coagulation between 20-25 or less near the nerve root 10-15. High temperature can be caused for high intensity power and is risky for tissue damage. Respecting the recommended power limits avoids potential

complications. In both systems a short-term and brief use is recommended rather than one prolonged application.

As a possibility, available in other centers, the radiofrequency system is also frequently used for control the hemorrhage during the surgery. The use of bone wax its also a possibility for specific situations. The coagulation devices used during UBE (Figure 4).



Figure 4: A: monopolar device; B: bipolar devices; C: intraoperative image with bipolar cautery.

During the surgery the fluid quality and type, temperature, infusion pressure and also the input and out put of the fluid are very important.

As described above the UBE its a surgery performed "underwater" and the correct circulation of the fluid in the surgical field permit a better vision and consequently a correct identification of anatomical structures and can avoid many complications. In our center we use a saline solution infused by gravity, with a infusion rate less than 30 mmHg.

All the procedures are under general anesthesia, we consider the benefit much more high then the risk. In other centers they use different methods such as local, epidural, and also general anesthesia.

The surgical procedures are performed in a specific operating table and support device for spinal surgery. The distribution in our operating room, the surgeon are located in the side of the mainly symptom presented by the patient and the endoscopic system in front of the surgeon. The irrigating saline lines in the distal portion of the surgical table.

In the case of a right-handed surgeon, the left hand the non-dominant holds the endoscope. In this case the endoscope-related lines, camera cable, and irrigating saline line, should be tucked to the left side of the operator. The drilling system should be fixed at the right side of the surgeon in left sides approaches.

After after anesthetic induction a very important point is the patient position, the first step of the surgical procedure. A correct prone position avoid

Auctores Publishing LLC – Volume 17(2)-358 www.auctoresonline.org ISSN: 2578-8868 complications and is mandatory check all the points and the correct distribution of the patient body, to obtain no pressure in eyes, and reduce in all articulations. The abdominal hight pressure will cause more bleeding, particularly in epidural space. Reducing the pressure to this level can make the entire procedure easier.

We recommend also flex the legs to indirectly open the vertebral space, very important particularly for interlaminal approaches.

The placement of surgical drapes, with all care to ensure asepsis, is equally important taking into account that the UBE surgery use many bags of saline fluid to irrigation. In this context ideally a kit with a collection bag included should be available. If not we created a simple aseptic working field of 10 cm around the skin incisions with the surgical drapes and place a antimicrobial film. After that we make a bulge around the surgical site with compresses and place sterilized bags on the side to collect the used saline solution. Finally we also recommend the application of a second antimicrobial film to isolate the surgical field.

Related with saline solution irrigation, the hydraulic pressure is a important point for a safe and smooth operation, in our case we can controlled this point using the gravity as described above. We placed the saline solution bag 50–70 cm above the patient's back, for obtain a goal of water pressure inside the body around 30 mmHg.

J. Neuroscience and Neurological Surgery

Copy rights @ Gonçalo Januário,

The UBE procedure use many instruments common to open spine surgery, is a real advantage for the hospitals because represent a reduction in the costs if we compare with uniportal endoscopy.

About the surgical instruments as commented previously in our center we use a common instrument but we also have exclusive instruments for UBE.

Actually we have also available many specially designed UBE toolset, which is comprised common tools and also devices for establishing the standard UBE surgical pathway and performing all surgical steps. Particular instruments as rotating Kerrison punch, curved Kerrison punch, curved Curets and the Indian knife all of them are essential for do UBE surgery (Figure 5).



Figure 5 : A: Kerrison punch ; B: curets; C: Indian knife, D: Pituitary rongeurs.

Many of the surgical advances that have emerged over the time are dependent and related to anesthetic advances. Like in all the surgical procedures the blood pressure is a important parameter, in UBE surgery is a crucial point and a optimal blood pressure (BP) is essential for a bloodless endoscopic view. For this endoscopic procedures we recommend a optimal systolic BP between 90–100 mmHg in general anesthesia, based in our experience. Further randomized and systemized study is needed.

For perform all surgical techniques and particularly for UBE is mandatory a complete understanding of the Basic Concepts. An adequate knowledge of the anatomy of the area to be operated is the main concept, and after all the learning curve like for all the surgical techniques.

At the beginning its mandatory understand and utilize the new basic concepts, take care of each of the steps in order to reduce the rate of complications, which is usually higher in this phase of the training process. These basic concepts should naturally include knowledge of the function of the different instruments used during the procedure. Every surgeon has their own learning curve, it is up to each one to do the best possible to provide greater benefit to the patient with each particular technique. By doing so, a consistently desirable surgical outcome is guaranteed while minimizing the possibility of complications.

The quality of the learning curve are related with many factors for examples the number of cases in a period of time, a professor for clarify any type of doubts, previous knowledges and personal experience.

The pathology to be treated, with indication for surgical treatment, will determine the use of one approach or another one, adjusting the technique to the patient's needs.

It is equally important to know the objectives to be achieved with the technique as well as its limitations, complications or controversies.

The majority of spine surgeons have experience with many different methods for treat the different types of spine pathologies, and they use the most appropriate for each specific cases. Frequently they often use the one they know better and have high control and fewer complications.

All surgical technique present advantages and disadvantages and has its own distinct basic concepts.

Basic Concepts of UBE:

1. Two-Port Approach:

- One working portal for inserting instruments (like probes, forceps, or drills).
- One viewing portal for a high-definition camera or endoscope that

Auctores Publishing LLC – Volume 17(2)-358 www.auctoresonline.org ISSN: 2578-8868

provides the surgeon with a clear visualization of the surgical site.

2. Unilateral Access:

• Both portals are on the same side of the patient, reducing the need for muscle dissection and retraction on the opposite side. This decreases soft tissue trauma.

3. Endoscopic Visualization:

- High-definition endoscopes with fiber-optic technology provide a clear, magnified view of the surgical area.
- 3D imaging and navigation systems may also be used to enhance precision.

4. Minimally Invasive:

 Smaller incisions (typically 1–2 cm each) are used, which leads to reduced blood loss, shorter recovery times, and less postoperative pain compared to open surgery.

5. Indications:

- Primarily used in spinal surgery (e.g., discectomy, laminectomy, and spinal fusion).
- It can also be applied in other areas, such as orthopedic or soft tissue procedures, when appropriate.

Advantages of UBE:

- Reduced Recovery Time: Faster postoperative recovery due to smaller incisions and less disruption of tissues.
- Less Blood Loss: The minimal invasiveness of the procedure results in less bleeding.
- Reduced Risk of Infection: Smaller incisions and less tissue exposure help minimize infection risk.
- Better Cosmetic Outcome: Smaller scars or sometimes no visible scarring depending on the procedure.

Challenges:

- Steeper Learning Curve: Requires significant training and skill to operate endoscopic equipment effectively.
- Limited Working Space: The confined workspace and specialized tools require precision and careful planning.

• Cost: The high-tech equipment involved can make the procedure more expensive than traditional approaches.

UBE is an evolving field, and with advancements in technology, its applications are becoming more widespread in spinal and other types of surgeries.

In resume we divide the UBE technique in 5 basic concepts as follows:

1. Unilateral biportal endoscopy

Previously called biportal endoscopic spine surgery (BESS) name related with the use of two portals for do the surgery. Actually the nomenclature in use is Unilateral biportal endoscopy (UBE).

This technique applies the knee arthroscopy surgical concept to the spine. In this approach we made two portals in the same side, unilateral, of the spine pathology. One of them is the endoscopic portal and the other is the working portal for introduce all the materials that we need during the procedure.

The greats benefits of the UBE are that uses the same surgical route as in open surgery, and also very important the movements of the endoscopic and the different tools have independent ways in separate portals. This point reduce visual limitations and permit a high range of movement for all the tools and devices.

2. Surgery in liquid medium

The UBE uses irrigating saline fluid. This solution in comparison with other increase the conduction and finally the high precision on the homeostasis devices.

During the surgical procedure, the work space is filled with saline solution that enters through the endoscope portal. Continuous irrigation is very important, is essential control the input and output of the solution to balance the hydrostatic pressure in the work space. Its important control the pressor inside the canal for reduce the complications rate.

This balance is critical for maintain a clear vision of the surgical field and facilitating the surgery and turn on more safety. It has also been proven that continuous fluid is associated with a lower rate of infections.

3. Working movements in Triangulation

In the endoscope the lens are located in the front position of the device. In the beginning stage of the learning curve it's hard to place the different instruments under endoscopic view.

During the procedure, the endoscope and the different instruments must be moved in a triangular or convergent manner. In this way the distal end of each instrument is visible allowing safe work. Triangulation means placing the ends of the instruments just in front of the endoscope lens. If the skin incision are made in the correct place the endoscope and the different tools don't cross or shock. When this situation occur and the endoscope and instruments are crossing each other in the shaft, each interferes with the other's movement complicating the performance of the procedure.

This situation is named 'early triangulation' and in this situation is mandatory or recommended a repositioning, using for them a c-arm or Oarm fluoroscopy or with neuronavigation.

The most extreme situation would be when the endoscope cannot face the instruments at the tip, it will not be possible to find the instrument in the endoscopic view. This situation is called "open triangulation". The opening of the fascia should be reviewed and, if necessary, repositioning should also be carried out, using the intraoperative imaging devices available in each operating room.

The UBE is a minimal invasive approach the muscle destruction is reduced, during the approach a important point is a correct open the muscular fascia. In our case we open on cruce, in both portals, and after Auctores Publishing LLC – Volume 17(2)-358 www.auctoresonline.org ISSN: 2578-8868

with a muscle detacher create a cavity through the soft tissues and multifidus muscle. This initial step is very important and guarantees better output flow and a clear vision.

With this approach we obtain a free movement of the different devices and tools with less destruction, and safety, with small skin incisions, less then 1 cm.

4. Two hands procedure, the coordination.

In UBE the non-dominant hand supports the main acting hand for obtain safe performances. In this technique the non-dominant hand holds the endoscope, and the dominant hand operates the surgical instruments. The dominant hand should function without any support from the nondominant hand. The most important point is do the surgery safely and respect the ergonomic recommendations.

To access certain points, blunt angled instruments or changing the angulation of the endoscope lens are sometimes necessary, as 30° .

Surgical skills must be developed to obtain the division and elevation technique, while removing the yellow ligament and ensuring sufficient surgical space, performing all movements smoothly to achieve a safe surgery particularly when we are inside the spinal canal, close to the spinal canal, near the duramater or the nerve roots.

In UBE the instrument and endoscope should be handled with pivot movements. This type of movements should be made, using the two portals as the central pivot point. As an example, if you plan to move the tip of the instruments forward, then you should move the tip of the other instruments, which are outside of the patient's body, backward. This theory applies to every movement performed during the procedure.

This is a specific knowledge that you need to introduce and become accustomed to this during the learning curve period.

5. Image quality

After introducing the endoscope inside the patient's body, the distance to the area to be treated is reduced, this is one of the great differences and advantages in relation to surgery performed with a microscope. The advantage is even more evident in cases of bilateral approaches in which by entering from one side, in cases of interlaminar approaches, it is possible to reach the contralateral side, achieving adequate central canal and lateral recess decompression with release of the nerve roots in the contralateral foramen.

For obtain a good vision its also important that the irrigation saline solution flows continuously through the tip of the endoscope sheath, and the flow pressure controls bleeding to some extent and removes hematoma to leave a clear surgical field.

For obtain a estable vision without trembling caused by the surgeon harm, did placing the arm close to the body.

Working in triangulation formation, as described in the point 3, the endoscope and the instruments work independently but very adjacent to each other this point permit a wide vision but exit the risk of causing damage to the endoscope lens particularly by speeding drill burrs. To avoid this situation its important maintain a distance between the endoscope and the instruments.

As recommendation remember that a small movement (rotation, introduction and extraction of the endoscope) of the tip of the endoscope results in very different surgical vision. The UBE technique is the present and the future for the treatment of vertebral pathologies, in the different spine levels. It's available in more and more places around the world and maybe soon will be the gold standard technique for the treatment of disc herniation and other degenerative diseases. Mainly in stand alone decompression surgeries and in other cases associate with minimal invasive, percutaneous fusion techniques.

J. Neuroscience and Neurological Surgery

Acknowledges:

All the patients, all my colleagues during my career especially to my mentor in endoscopic spine surgery Dr. Pedro Llinás. To my family for all the support.

References:

- Spaner SJ, Warnock GL. A brief history of endoscopy, laparoscopy, and laparoscopic surgery. J Laparoendosc Adv Surg Tech A. 1997;7(6):369–73.
- 2. Yong Ahn, Semin Lee; Uniportal versus biportal endoscopic spine surgery: a comprehensive review, Volume 20, 2023 -

Issue 7, Pages 549-556 | Received 25 Jul 2022, Accepted 12 May 2023, Published online: 15 May 2023.

- 3. Mayer HM. A history of endoscopic lumbar spine surgery: what have we learnt? Biomed Res Int. 2019;2019:4583943.
- 4. Sang-Kyu Son, Dong Han Kim, and Hayati Aygün; The Basics and Concepts of Unilateral Biportal Endoscopy; Springer Nature Singapore Pte Ltd. 2022 9 D. H. Heo et al. (eds.), Unilateral Biportal Endoscopic Spine Surgery.
- Arangio G, Kostelnik KE. Intraarticular pressures in a gravityfed arthroscopy fluid delivery system. Arthroscopy. 1992;8(3):341–4.



This work is licensed under Creative Commons Attribution 4.0 License

To Submit Your Article Click Here:

Submit Manuscript

DOI:10.31579/2578-8868/358

Ready to submit your research? Choose Auctores and benefit from:

- ➢ fast, convenient online submission
- > rigorous peer review by experienced research in your field
- rapid publication on acceptance
- > authors retain copyrights
- > unique DOI for all articles
- immediate, unrestricted online access

At Auctores, research is always in progress.

Learn more https://auctoresonline.org/journals/neuroscience-and-neurologicalsurgery